

Nur für Forschungszwecke

Neurocan Polyklonaler Antikörper

Katalog-Nr.: 27087-1-AP



Allgemeine Informationen

Katalog-Nr.:
27087-1-AP

Größe:

150ul, Konzentration: 1000 µg/ml von 1463

Nanodrop und 567 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG25821

GenBank-Zugangsnummer:

NM_004386

GeneID (NCBI):

1463

Vollständiger Name:

neurocan

Berechnete Masse:

143 kDa

Beobachtete Masse:

270 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:1000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human

Positivkontrollen:

WB : fetales humanes Hirngewebe,

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Fetal human brain lysates were subjected to SDS PAGE followed by western blot with 27087-1-AP (Neurocan antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.