

NKCC1,SLC12A2 Polyklonaler Antikörper

Katalog-Nr.:28055-1-AP


Allgemeine Informationen

Katalog-Nr.: 28055-1-AP	GenBank-Zugangsnummer: BC033003	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul , Konzentration: 1000 µg/ml von6558 Nanodrop;	GeneID (NCBI): 15058	Empfohlene Verdünnungen: WB 1:2000-1:10000 IHC 1:300-1:1200 IF 1:50-1:500
Wirt: Kaninchen	Vollständiger Name: solute carrier family 12 (sodium/potassium/chloride transporters), member 2	
Isotyp: IgG	Berechnete Masse: 131 kDa	
Immunogen Katalognummer: AG27895	Beobachtete Masse: 160 kDa	


Anwendungen

Geprüfte Anwendungen: IF, IHC, WB, ELISA	Positivkontrollen:
Getestete Reaktivität: Human	WB : ungekochte HeLa-Zellen,
Hinweis-IHC: Antigendemaskierung mit TE-Puffer pH 9,0 empfohlen	


Ausgewählte Validierungsdaten



HeLa cells were subjected to SDS PAGE followed by western blot with 28055-1-AP (SLC12A2 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human breast cancer tissue slide using 28055-1-AP (SLC12A2 antibody) at dilution of 1:600 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human breast cancer tissue slide using 28055-1-AP (SLC12A2 antibody) at dilution of 1:600 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).

