

Nur für Forschungszwecke

FC tag Polyklonaler Antikörper

Katalog-Nr.: 28108-1-AP



Allgemeine Informationen

Katalog-Nr.:

28108-1-AP

Größe:

150ul, Konzentration: 150 µg/ml von
Nanodrop;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG27844

GenBank-Zugangsnummer:

GeneID (NCBI):

Vollständiger Name:

Berechnete Masse:

26 kDa

Beobachtete Masse:

55 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:2000-1:10000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human

Positivkontrollen:

WB : Recombinant protein,

Hintergrundinformationen

This antibody detects the Fc fragment of mouse IgG2a.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 28108-1-AP (FC tag antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.