

TRIM3 Polyklonaler Antikörper

Katalog-Nr.: 28392-1-AP

Allgemeine Informationen

Katalog-Nr.:

28392-1-AP

Größe:

150ul, Konzentration: 700 µg/ml von Nanodrop und 333 µg/ml durch die Bradford-Methode mit BSA als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG28823

GenBank-Zugangsnummer:

NM_001248006

GeneID (NCBI):

10612

Vollständiger Name:

tripartite motif-containing 3

Berechnete Masse:

81 kDa

Beobachtete Masse:

75 kDa, 80-85 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:1000-1:6000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Positivkontrollen:

WB: Maushirngewebe, Rattenhirngewebe

Hintergrundinformationen

TRIM3, a member of the tripartite motif (TRIM) family, also called the 'RING-B-box-coiled-coil' (RBCC) subgroup of RING finger proteins. The TRIM family of genes is largely studied because of their roles in development, differentiation and host cell antiviral defenses. TRIM3 is a mammalian tumor suppressor whose loss can increase the incidence and promote the development of GBM (PMID: 23318451). TRIM3 has four isoforms with the MW of 67-80 kDa and can be detected as 75, 80, 82, 85 kDa after posttranslational modification (PMID: 24947043/PMID: 24086586).

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 28392-1-AP (TRIM3 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.