

Allgemeine Informationen

Katalog-Nr.: 29440-1-AP	GenBank-Zugangsnummer: BC146793	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul, Konzentration: 600 µg/ml von Nanodrop;	GeneID (NCBI): 5364	Empfohlene Verdünnungen: WB 1:500-1:2000 IHC 1:50-1:500
Wirt: Kaninchen	Vollständiger Name: plexin B1	
Isotyp: IgG	Berechnete Masse: 2135 aa, 232 kDa	
Immunogen Katalognummer: AG31224	Beobachtete Masse: 200-250 kDa	

Anwendungen

Geprüfte Anwendungen: IHC, WB, ELISA	Positivkontrollen:
Getestete Reaktivität: Human, Maus, Ratte	WB : Maushirngewebe, Rattenhirngewebe IHC : Maushirngewebe,

Hinweis-IHC: Antigendemaskierung mit TE-Puffer pH 9,0 empfohlen. (*) Wahlweise kann die Antigendemaskierung auch mit Citratpuffer pH 6,0 erfolgen.

Hintergrundinformationen


Semaphorins, originally identified as axon guidance molecules, have also been implicated in angiogenesis, immunoregulation, and cancer (PMID: 16584533). Plexins are main receptors for semaphorins and are thought to control many of the functional effects of semaphorins. PLXNB1 is one of the three plexin-B family members. It is a high-affinity receptor for semaphorin 4D (SEMA4D). Interaction of SEMA4D and PLXNB1 has important physiologic effects in the immune and nervous systems, and is also involved in tumor progression, especially in tumor angiogenesis (PMID: 20858260).

Lagerung


Lagerungsbedingungen:
Bei -20°C lagern.
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 29440-1-AP (PLXNB1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using 29440-1-AP (PLXNB1 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).