

Allgemeine Informationen

Katalog-Nr.:
29751-1-AP

Größe:
150ul, Konzentration: 450 µg/ml von
Nanodrop;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG31335

GenBank-Zugangsnummer:
BC007545

GeneID (NCBI):
10273

Vollständiger Name:
STIP1 homology and U-box
containing protein 1

Beobachtete Masse:
35 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:
WB 1:1000-1:6000

Anwendungen

Geprüfte Anwendungen:
WB, ELISA

Getestete Reaktivität:
Human, Maus, Ratte

Positivkontrollen:

WB: A549-Zellen, HEK-293-Zellen, HeLa-Zellen,
Mausnierengewebe, NIH/3T3-Zellen, PC-12-Zellen,
Rattennierengewebe

Hintergrundinformationen

STIP1 Homology and U-Box Containing Protein 1 (STUB1), formerly known as CHIP (C-terminus of HSP70-Interacting Protein), is a chaperone-dependent E3 ubiquitin ligase. STUB1 is a 34.5 kDa co-chaperone encoded by the STUB1 gene. STUB1 is highly conserved across species and is expressed in most tissues.

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern.

Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 29751-1-AP (STUB1 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.