

Nur für Forschungszwecke

DOCK5 Polyklonaler Antikörper

Katalog-Nr.: 29938-1-AP



Allgemeine Informationen

Katalog-Nr.:

29938-1-AP

Größe:

150ul, Konzentration: 700 µg/ml von

Nanodrop;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG32193

GenBank-Zugangsnummer:

BC137175

GeneID (NCBI):

80005

Vollständiger Name:

dedicator of cytokinesis 5

Berechnete Masse:

1870 aa, 215 kDa

Beobachtete Masse:

190 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:1000-1:4000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Positivkontrollen:

WB: Jurkat-Zellen, NIH3T3-Zellen, Rattenhirngewebe

Hintergrundinformationen

DOCK5 is a member of the dedicator of cytokinesis protein family. Members of this family act as guanine nucleotide exchange factors for small Rho family G proteins. DOCK5 is thought to associate with adaptors CRK and CRKL, and function in regulation of intestinal epithelial cell spreading and migration on collagen IV.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 29938-1-AP (DOCK5 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.