

Nur für Forschungszwecke

INPP5J Polyklonaler Antikörper

Katalog-Nr.: 29949-1-AP



Allgemeine Informationen

Katalog-Nr.:

29949-1-AP

Größe:

150ul, Konzentration: 450 µg/ml von

Nanodrop;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG29242

GenBank-Zugangsnummer:

BC109288

GeneID (NCBI):

27124

Vollständiger Name:

inositol polyphosphate-5-phosphatase J

Berechnete Masse:

1006 aa, 107 kDa

Beobachtete Masse:

70 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:3000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human, Maus

Positivkontrollen:

WB: A549-Zellen, HeLa-Zellen, Mausherzgewebe, MCF-7-Zellen

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 29949-1-AP (INPP5) antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.