

Nur für Forschungszwecke

HLA-DQA2 Polyklonaler Antikörper

Katalog-Nr.: 55072-1-AP



Allgemeine Informationen

Katalog-Nr.:
55072-1-AP

Größe:

150ul, Konzentration: 1000 µg/ml von 3118

Nanodrop und 507 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:

Kaninchen

Isotyp:

IgG

GenBank-Zugangsnummer:

NM_020056

GeneID (NCBI):

3118

Vollständiger Name:

major histocompatibility complex,
class II, DQ alpha 2

Berechnete Masse:

28 kDa

Beobachtete Masse:

28 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:200-1:1000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human, Maus

Positivkontrollen:

WB: Raji-Zellen, Mausmilzgewebe, Maus-
Thymusgewebe

Hintergrundinformationen

HLA-DQA2, also named as HLA-DXA, belongs to the MHC class II family. It binds peptides derived from antigens that access the endocytic route of antigen presenting cells (APC) and presents them on the cell surface for recognition by the CD4 T-cells. The antibody is specific to HLA-DQA2.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern.

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Ausgewählte Validierungsdaten



Raji cells were subjected to SDS PAGE followed by western blot with 55072-1-AP (HLA-DQA2 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.