

Nur für Forschungszwecke

Ly-6G Monoklonaler Antikörper

Katalog-Nr.: 65078-1-Ig



Allgemeine Informationen

Katalog-Nr.:

65078-1-Ig

Größe:

100 µg, 500 µg/ml

Wirt:

Ratte

Isotyp:

IgG2a, kappa

GenBank-Zugangsnummer:

X70920

GeneID (NCBI):

546644

Vollständiger Name:

Lymphocyte antigen 6 complex, locus G

Reinigungsmethode:

Affinitätsreinigung

CloneNo.:

1A8

Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Maus

Hintergrundinformationen

Ly-6G (lymphocyte antigen 6 complex, locus G), also known as Gr-1, is a 21-25 kDa, glycosylphosphatidylinositol-anchored protein expressed on myeloid lineage cells in mouse bone marrow (PMID: 8360469). The expression of Ly-6G increases on neutrophils as they differentiate from immature cells in the bone marrow to mature cells in the blood and spleen (PMID: 8890901). Antibodies targeting Ly6G (RB6-8C5 or 1A8) are commonly used in studies aimed at identifying the role of neutrophils (PMID: 23543767). The 1A8 mAb is specific for Ly-6G (PMID: 8360469).

Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Stable for one year after shipment.

Lagerungspuffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ mouse bone marrow cells were surface stained with 0.5 ug Anti-Mouse Ly-6G (65078-1-Ig, Clone:1A8) and FITC anti-Rat IgG2a Antibody at dilution 1:100 (red), or 0.5 ug Control Antibody. Cells were not fixed.



1X10⁶ mouse bone marrow cells were surface stained with 0.5 ug Anti-Mouse Ly-6G (65078-1-Ig, Clone:1A8) and FITC anti-Rat IgG2a Antibody at dilution 1:100, or 0.5 ug Control Antibody. Cells were not fixed.