

Nur für Forschungszwecke

STRO-1 Monoklonaler Antikörper

Katalog-Nr.: 65184-1-Ig



Allgemeine Informationen

Katalog-Nr.:

65184-1-Ig

Größe:

100ug, 0.5 mg/ml

Wirt:

Maus

Isotyp:

IgM, lambda

GenBank-Zugangsnummer:

GeneID (NCBI):

Vollständiger Name:

Reinigungsmethode:

Mit Präzipitation und
Chromatographie aufgereinigt

CloneNo.:

STRO-1

Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Human

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Stable for one year after shipment.

Lagerungspuffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ lysed human blood cells were surface stained with 0.2 ug Anti-Human STRO-1 (65184-1-Ig, Clone:STRO-1) and AlexaFluor-Goat anti-Mouse IgG at dilution 1:1000. Cells were not fixed. CD45 dim, CD34 positive cells were gated.



1X10⁶ lysed human blood cells were surface stained with 0.2 ug Anti-Human STRO-1 (65184-1-Ig, Clone:STRO-1) and AlexaFluor-Goat anti-Mouse IgG at dilution 1:1000. Cells were not fixed. CD45 dim, CD34 positive cells were gated.