

Nur für Forschungszwecke

CD45RA Monoklonaler Antikörper

Katalog-Nr.: 65234-1-Ig



Allgemeine Informationen

Katalog-Nr.:

65234-1-Ig

Größe:

100ug, 0.5 mg/ml

Wirt:

Maus

Isotyp:

IgG1, kappa

GenBank-Zugangsnummer:

XM_006249910

GeneID (NCBI):

24699

Vollständiger Name:

protein tyrosine phosphatase,
receptor type, C

Berechnete Masse:

143 kDa

Reinigungsmethode:

Affinitätsreinigung

CloneNo.:

OX-33

Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Ratte

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Stable for one year after shipment.

Lagerungspuffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ Wistar rat splenocytes were surface stained with 0.5 μ g Anti-Rat CD45RA (65234-1-Ig, Clone:OX-33) and FITC anti-Mouse IgG1 Antibody at dilution 1:100, or 0.5 μ g Control Antibody. Cells were not fixed.