

Nur für Forschungszwecke

CD94 Monoklonaler Antikörper

Katalog-Nr.: 65240-1-Ig



Allgemeine Informationen

Katalog-Nr.:

65240-1-Ig

Größe:

100ug, 0.5 mg/mL

Wirt:

Ratte

Isotyp:

IgG2a, kappa

GenBank-Zugangsnummer:

BC120854

GeneID (NCBI):

16643

Vollständiger Name:

killer cell lectin-like receptor,
subfamily D, member 1

Reinigungsmethode:

Affinitätsreinigung

CloneNo.:

18D3

Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Maus

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Stable for one year after shipment.

Lagerungspuffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ mouse splenocytes were surface stained with 1 μ g Anti-Mouse CD94 (65240-1-Ig, Clone:18D3) and FITC anti-rat IgG2a Antibody at dilution 1:100 (red), or 1 μ g Control Antibody. Cells were not fixed.