

Nur für Forschungszwecke

Nectin-4/PVRL4 Monoklonaler Antikörper

Katalog-Nr.: 67721-1-Ig



Allgemeine Informationen

Katalog-Nr.:

67721-1-Ig

Größe:

150ul, Konzentration: 1000 µg/ml von 81607

Nanodrop und 466 µg/ml durch die Bradford-Methode mit BSA als Standard;

Wirt:

Maus

Isotyp:

IgG1

Immunogen Katalognummer:

AG14414

GenBank-Zugangsnummer:

BC010423

GeneID (NCBI):

81607

Vollständiger Name:

poliovirus receptor-related 4

Berechnete Masse:

510 aa, 55 kDa

Beobachtete Masse:

43 kDa

Reinigungsmethode:

Protein-G-Reinigung

CloneNo.:

2D4B6

Empfohlene Verdünnungen:

WB 1:500-1:3000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human

Positivkontrollen:

WB : humanes Plazenta-Gewebe, A431-Zellen

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern.

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



human placenta tissue were subjected to SDS PAGE followed by western blot with 67721-1-Ig (Nectin-4/PVRL4 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.