

Nur für Forschungszwecke

PSIP1 Monoklonaler Antikörper

Katalog-Nr.:67723-1-Ig



Allgemeine Informationen

Katalog-Nr.: 67723-1-Ig	GenBank-Zugangsnummer: NM_021144	Reinigungsmethode: Protein-G-Reinigung
Größe: 150ul	GeneID (NCBI): 11168	CloneNo.: 1E9E7
Wirt: Maus	Vollständiger Name: PC4 and SFRS1 interacting protein 1	
Isotyp: IgG1	Berechnete Masse: 60 kDa	
Immunogen Katalognummer: AG30562		

Anwendungen

Geprüfte Anwendungen:
ELISA
Getestete Reaktivität:
Human, Maus

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern.
Lagerungspuffer:
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

