

Nur für Forschungszwecke

Human IgG Kappa chain Monoklonaler Antikörper

Katalog-Nr.: 67761-1-Ig



Allgemeine Informationen

Katalog-Nr.:

67761-1-Ig

Größe:

150ul, Konzentration: 1000 µg/ml von
Nanodrop und 439 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:

Maus

Isotyp:

IgG1

GenBank-Zugangsnummer:

GeneID (NCBI):

Vollständiger Name:

Berechnete Masse:

155 kDa

Beobachtete Masse:

~25 kDa

Reinigungsmethode:

Protein-A-Reinigung

CloneNo.:

4D11F7

Empfohlene Verdünnungen:

WB 1:5000-1:50000

IHC 1:2500-1:10000

Anwendungen

Geprüfte Anwendungen:

IHC, WB, ELISA

Getestete Reaktivität:

Human

Positivkontrollen:

WB : human plasma,

IHC : humanes Tonsillitisgewebe,

**Hinweis-IHC: Antigendemaskierung mit TE-
Puffer pH 9,0 empfohlen. (*) Wahlweise
kann die Antigendemaskierung auch mit
Citratpuffer pH 6,0 erfolgen.**

Hintergrundinformationen

This antibody detects the Kappa chain of human IgG,

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern.

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.


Ausgewählte Validierungsdaten



Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 67761-1-Ig (Human IgG Kappa chain antibody) at dilution of 1:5000 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 67761-1-Ig (Human IgG Kappa chain antibody) at dilution of 1:5000 (under 4x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 67761-1-Ig (Human IgG Kappa chain antibody) at dilution of 1:5000 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



16 nL human plasma were subjected to SDS PAGE followed by western blot with 67761-1-Ig (Human IgG Kappa chain antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.