

Nur für Forschungszwecke

NSF Monoklonaler Antikörper

Katalog-Nr.: 68061-1-Ig

Vorgestelltes Produkt



Allgemeine Informationen

Katalog-Nr.:

68061-1-Ig

Größe:

150ul, Konzentration: 1000 µg/ml von 4905

Nanodrop;

Wirt:

Maus

Isotyp:

IgG2a

Immunogen Katalognummer:

AG15422

GenBank-Zugangsnummer:

BC030613

GeneID (NCBI):

von 4905

Vollständiger Name:

N-ethylmaleimide-sensitive factor

Berechnete Masse:

744 aa, 83 kDa

Beobachtete Masse:

70-83 kDa

Reinigungsmethode:

Protein-A-Reinigung

CloneNo.:

2B6H6

Empfohlene Verdünnungen:

WB 1:5000-1:50000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Hausschwein, Huhn, Human, Kaninchen, Maus, Ratte

Positivkontrollen:

WB: Hausschwein-Hirngewebe, HEK-293-Zellen, Hühnerhirngewebe, Kaninchenhirngewebe, Maushirngewebe, PC-12-Zellen, Rattenhirngewebe, SH-SY5Y-Zellen

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Ausgewählte Validierungsdaten




Various lysates were subjected to SDS PAGE followed by western blot with 68061-1-Ig (NSF antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



PC-12 cells were subjected to SDS PAGE followed by western blot with 68061-1-Ig (NSF antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



WB result of NSF antibody (68061-1-Ig; 1:60000; incubated at room temperature for 1.5 hours) with sh-Control and sh-NSF transfected HEK-293 cells.



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with 68061-1-Ig (NSF antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.