

Allgemeine Informationen

Katalog-Nr.:

81290-1-RR

GenBank-Zugangsnummer:

GeneID (NCBI):

Vollständiger Name:

Reinigungsmethode:

Protein-A-Reinigung

CloneNo.:

6I21

Empfohlene Verdünnungen:

WB 1:5000-1:50000

IF 1:375-1:1500

Größe:

100ul, Konzentration: 750 µg/ml von

Nanodrop;

Wirt:

Kaninchen

Isotyp:

IgG

Anwendungen

Geprüfte Anwendungen:

IF, WB, ELISA

Getestete Reaktivität:

Rekombinantes Protein

Positivkontrollen:

WB : Transfizierte HEK-293-Zellen,

IF : Transfizierte HEK-293-Zellen,

Hintergrundinformationen

Protein tags are a protein or peptide sequences located either on the C- or N- terminal of the target protein, which can facilitate solubility, detection, purification, localization, and expression of the target protein. The HA tag corresponds to amino acid residues YPYDVPDYA of a surface glycoprotein -human influenza virus hemagglutinin (HA). The HA tag is commonly used for a variety of research applications including chromatin immunoprecipitation, ELISA, flow cytometry, western blotting, immunocytochemistry/immunofluorescence among of others.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Ausgewählte Validierungsdaten



HA-tag recombinant protein, non-transfected HEK-293 cells, transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with 81290-1-RR (HA Tag antibody) at dilution of 1:30000 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (-20°C Ethanol) fixed Transfected HEK-293 cells using HA Tag antibody (81290-1-RR, Clone: 6I21) at dilution of 1:750 and CoraLite®594-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



IP result of anti-HA Tag(IP:81290-1-RR, 4ug; Detection:81290-1-RR 1:20000) with Transfected HEK-293T cells lysate 795 ug.