

Nur für Forschungszwecke

STRO-1 Monoklonaler Antikörper

Katalog-Nr.: AB-65184



Allgemeine Informationen

Katalog-Nr.:

AB-65184

Größe:

100tests, 5 µl/test

Wirt:

Maus

Isotyp:

IgM, lambda

GenBank-Zugangsnummer:

GeneID (NCBI):

Vollständiger Name:

Reinigungsmethode:

Mit Präzipitation und
Chromatographie aufgereinigt

CloneNo.:

STRO-1

Anregungs-/Emissionsmaxima-
Wellenlängen:

404 nm / 458 nm

Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Human

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:

PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ human whole blood cells were surface stained with 5 ul Atlantic Blue™ Anti-Human STRO-1 (AB-65184, Clone: STRO-1) or unstained. Cells were not fixed. CD45dim, CD34+ cells were gated.