

Nur für Forschungszwecke

IgG2a Isotype Control Monoklonaler Antikörper



Katalog-Nr.:CL405-65208

Allgemeine Informationen

Katalog-Nr.:
CL405-65208
Größe:
100ug, 0.5 mg/ml
Wirt:
Maus
Isotyp:
IgG2a, kappa

GenBank-Zugangsnummer:
GeneID (NCBI):
Vollständiger Name:

Reinigungsmethode:
Affinitätsreinigung
CloneNo.:
C1.18.4
Anregungs-/Emissionsmaxima-Wellenlängen:
399 nm / 422 nm

Anwendungen

Geprüfte Anwendungen:
FC
Getestete Reaktivität:

Hintergrundinformationen

The C1.18.4 immunoglobulin is useful as an isotype-matched control. The C.1.18.4 immunoglobulin has an unknown binding specificity and is used as an isotype control for mouse IgG2a antibodies. This antibody has been quality-tested for flow cytometry as negative control.

Lagerung

Lagerungsbedingungen:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Lagerungspuffer:
PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ mouse splenocytes were surface stained with 0.5 ug CoraLite® Plus 405 Mouse IgG2a Isotype Control (CL405-65208, Clone: C1.18.4) (red), or 0.5 ug Control Antibody (green), or unstained (blue). Cells were not fixed.