

Nur für Forschungszwecke

CD94 Monoklonaler Antikörper

Katalog-Nr.: CL405-65240



Allgemeine Informationen

Katalog-Nr.:

CL405-65240

Größe:

100ug , 500 µg/ml

Wirt:

Ratte

Isotyp:

IgG2a, kappa

GenBank-Zugangsnummer:

BC120854

GeneID (NCBI):

16643

Vollständiger Name:

killer cell lectin-like receptor,
subfamily D, member 1

Reinigungsmethode:

Affinitätsreinigung

CloneNo.:

18D3

Anregungs-/Emissionsmaxima-
Wellenlängen:

399 nm / 422 nm

Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Maus

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ mouse splenocytes were surface stained with 0.5 ug CoraLite® Plus 405 Anti-Mouse CD94 (CL405-65240, Clone: 18D3) or 0.5 ug Control Antibody. Cells were not fixed.