

Nur für Forschungszwecke

NONO Polyklonaler Antikörper

Katalog-Nr.:CL488-11058

Vorgestelltes Produkt



Allgemeine Informationen

Katalog-Nr.:
CL488-11058

Größe:

100ul , Konzentration: 1000 µg/ml von4841

Nanodrop;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG1512

GenBank-Zugangsnummer:

BC010049

GeneID (NCBI):

von4841

Vollständiger Name:

non-POU domain containing,
octamer-binding

Berechnete Masse:

54 kDa

Beobachtete Masse:

54 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anregungs-/Emissionsmaxima-

Wellenlängen:

493 nm / 522 nm

Anwendungen

Geprüfte Anwendungen:

FC (Intra)

Getestete Reaktivität:

Human, Maus, Ratte

Hintergrundinformationen

Non-pou domain-containing octamer-binding protein(NONO), other named p54nrb, belongs to DBHS family, whose protein are predominately nuclear and involves in transcriptional and posttranscriptional gene regulatory functions as well as DNA repair. This family has roles in a wide range of biological processes, such as the regulation of circadian rhythm, carcinogenesis and progression of cancer. NONO forms heterodimer with SFPO, and has been implicated in nuclear retention of A- to I-edited RNA as paraspeckle components, pre-mRAN 3'-end formation, cAMP cycling and transcriptional activation.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Vor Licht schützen.

Lagerungspuffer:

BS mit 50% Glycerin, 0,05% Proclin300, 0,5% BSA, pH 7,3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1×10^6 HepG2 cells were intracellularly stained with 0.4 μg CoraLite® Plus 488 Anti-Human NONO (CL488-11058) (red), or 0.4 μg Control Antibody. Cells were fixed and permeabilized with Transcription Factor Staining Buffer Kit (PF00011).