

Nur für Forschungszwecke

STRO-1 Monoklonaler Antikörper

Katalog-Nr.:CL488-65184



Allgemeine Informationen

Katalog-Nr.:
CL488-65184
Größe:
100tests , 5 µl/test
Wirt:
Maus
Isotyp:
IgM, lambda

GenBank-Zugangsnummer:
GeneID (NCBI):
Vollständiger Name:

Reinigungsmethode:
Mit Präzipitation und
Chromatographie aufgereinigt
CloneNo.:
STRO-1
Anregungs-/Emissionsmaxima-
Wellenlängen:
493 nm / 522 nm

Anwendungen

Geprüfte Anwendungen:
FC
Getestete Reaktivität:
Human

Hintergrundinformationen


Lagerung

Lagerungsbedingungen:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Lagerungspuffer:
PBS with 0.09% sodium azide and 0.5% BSA.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ human peripheral blood platelets were surface stained with 5 μ l CoraLite® Plus 488 Anti-Human STRO-1 (CL488-65184, Clone:STRO-1) or Mouse IgM Isotype Control. Cells were not fixed. CD45 dim, CD34⁺ lymphocytes were gated.