

Nur für Forschungszwecke

# CD49b Monoklonaler Antikörper

Katalog-Nr.:CL488-65252



## Allgemeine Informationen

<b>Katalog-Nr.:</b> CL488-65252	<b>GenBank-Zugangsnummer:</b> NM_002203	<b>Reinigungsmethode:</b> Affinitätsreinigung
<b>Größe:</b> 100tests , 5 µl/test	<b>GeneID (NCBI):</b> 3673	<b>CloneNo.:</b> AK7
<b>Wirt:</b> Maus	<b>Vollständiger Name:</b> integrin, alpha 2 (CD49B, alpha 2 subunit of VLA-2 receptor)	<b>Anregungs-/Emissionsmaxima-Wellenlängen:</b> 493 nm / 522 nm
<b>Isotyp:</b> IgG1, kappa	<b>Berechnete Masse:</b> 126 kDa	

## Anwendungen

**Geprüfte Anwendungen:**  
FC

**Getestete Reaktivität:**  
Human

## Hintergrundinformationen

### Lagerung

**Lagerungsbedingungen:**  
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

**Lagerungspuffer:**  
PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Ausgewählte Validierungsdaten



1X10<sup>6</sup> human PBMCs were surface stained with 5 ul CoraLite® Plus 488 Anti-Human CD49b (CL488-65252, Clone:AK7) (red) or unstained. Cells were not fixed. Monocytes were gated.