

Nur für Forschungszwecke

FCGR2A / CD32a Monoklonaler Antikörper

Katalog-Nr.:CL488-66529



Allgemeine Informationen

Katalog-Nr.: CL488-66529	GenBank-Zugangsnummer: BC020823	Reinigungsmethode: Protein-A-Reinigung
Größe: 100ul , Konzentration: 1000 µg/ml von2212 Nanodrop;	GeneID (NCBI): Vollständiger Name:	CloneNo.: 1D3F1
Wirt: Maus	Berechnete Masse: 35 kDa	Empfohlene Verdünnungen: IF 1:50-1:500
Isotyp: IgG2b	Immunogen Katalognummer: AG25878	Anregungs-/Emissionsmaxima-Wellenlängen: 493 nm / 522 nm

Anwendungen

Geprüfte Anwendungen: FC (Intra), IF	Positivkontrollen: IF : humanes Tonsillitisgewebe,
Getestete Reaktivität: Human	

Hintergrundinformationen

IgG Fc receptors (FcγRs) play important roles in immune responses. Low affinity IgG Fc receptor Fc gamma RIla (FCGR2A, CD32a) is a cell surface receptor found on phagocytic cells such as macrophages and neutrophils, and is involved in the process of phagocytosis and clearing of immune complexes. This antibody is conjugated with CL488, Ex/Em 488 nm/515 nm.

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Vor Licht schützen. Nach dem Versand ein Jahr stabil.
Lagerungspuffer:
BS mit 50% Glycerin, 0,05% Proclin300, 0,5% BSA, pH 7,3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Ausgewählte Validierungsdaten



1x10⁶ K526 cells were stained with 0.20ug FCGR2A / CD32a antibody (CL488-66529, red) and control antibody (blue). Fixed with 90% MeOH.



Immunofluorescent analysis of (4% PFA) fixed human tonsillitis tissue using CL488-66529 (FCGR2A / CD32a antibody) at dilution of 1:50.



Immunofluorescent analysis of (4% PFA) fixed human tonsillitis tissue using CL488-66529 (FCGR2A / CD32a antibody) at dilution of 1:50.