

Nur für Forschungszwecke

Cyclin H Monoklonaler Antikörper

Katalog-Nr.:CL488-67065



Allgemeine Informationen

Katalog-Nr.: CL488-67065	GenBank-Zugangsnummer: BC005280	Reinigungsmethode: Protein-A-Reinigung
Größe: 100ul , Konzentration: 1000 µg/ml von902 Nanodrop;	GeneID (NCBI): 2C5G11	CloneNo.: 2C5G11
Wirt: Maus	Vollständiger Name: cyclin H	Empfohlene Verdünnungen: IF 1:50-1:500
Isotyp: IgG2a	Berechnete Masse: 38 kDa	Anregungs-/Emissionsmaxima- Wellenlängen: 493 nm / 522 nm
Immunogen Katalognummer: AG28556	Beobachtete Masse: 38 kDa	

Anwendungen

Geprüfte Anwendungen: IF	Positivkontrollen: IF : PC-3-Zellen,
Getestete Reaktivität: Human, Maus, Ratte	

Hintergrundinformationen

Lagerung


Lagerungsbedingungen:
Bei -20°C lagern. Vor Licht schützen.
Lagerungspuffer:
BS mit 50% Glycerin, 0,05% Proclin300, 0,5% BSA, pH 7,3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Immunofluorescent analysis of (4% PFA) fixed PC-3 cells using CoraLite® Plus 488 Cyclin H antibody (CL488-67065, Clone: 2C5G11) at dilution of 1:200, CL594-Phalloidin (red).