

SNX3 Monoklonaler Antikörper

Katalog-Nr.:CL488-68292

Allgemeine Informationen

Katalog-Nr.: CL488-68292	GenBank-Zugangsnummer: BC008444	Reinigungsmethode: Protein-G-Reinigung
Größe: 100ul , Konzentration: 1000 µg/ml von8724 Nanodrop;	GeneID (NCBI): von8724	CloneNo.: 2B5B7
Wirt: Maus	Vollständiger Name: sorting nexin 3	Empfohlene Verdünnungen: IF 1:50-1:500
Isotyp: IgG1	Berechnete Masse: 19 kDa	Anregungs-/Emissionsmaxima- Wellenlängen: 493 nm / 522 nm
Immunogen Katalognummer: AG1160	Beobachtete Masse: 16-20 kDa	

Anwendungen

Geprüfte Anwendungen: IF	Positivkontrollen: IF : A431-Zellen,
Getestete Reaktivität: Human, Maus, Ratte	

Hintergrundinformationen

SNX3 is a member of the sorting nexin family of proteins regulating membrane traffic. SNX3 regulates endosomal function through its PX (phox homology)-domain-mediated interaction with PtdIns(3)P (PMID: 11433298). It is a potential regulator of neurite formation and may play an important role in the development and function of the central nervous system (PMID: 19576982; 20817026).

Lagerung


Lagerungsbedingungen:
Bei -20°C lagern. Vor Licht schützen. Nach dem Versand ein Jahr stabil.
Lagerungspuffer:
BS mit 50% Glycerin, 0,05% Proclin300, 0,5% BSA, pH 7,3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Immunofluorescent analysis of (-20°C Methanol) fixed A431 cells using CoraLite® Plus 488 SNX3 antibody (CL488-68292, Clone: 2B5B7) at dilution of 1:200.