

Nur für Forschungszwecke

IgG2a Isotype Control Monoklonaler Antikörper



Katalog-Nr.:CL555-65244

Allgemeine Informationen

Katalog-Nr.:
CL555-65244
Größe:
100ug , 0.2 mg/mL
Wirt:
Maus
Isotyp:
IgG2a, kappa

GenBank-Zugangsnummer:
GeneID (NCBI):
Vollständiger Name:

Reinigungsmethode:
N/A
CloneNo.:
MOPC-173
Anregungs-/Emissionsmaxima-
Wellenlängen:
554 nm / 570 nm

Anwendungen

Geprüfte Anwendungen:
FC
Getestete Reaktivität:

Hintergrundinformationen


Lagerung

Lagerungsbedingungen:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Lagerungspuffer:
PBS with 0.09% sodium azide and 0.5% BSA.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ human PBMCs were surface stained with 0.2 ug CoraLite® Plus 555 Anti-Mouse IgG2a Isotype Control (CL555-65244, Clone:MOPC-173) (red), CL555 Anti-Human CD45RO (green), or unstained (blue). Cells were not fixed. Lymphocytes were gated.