

Nur für Forschungszwecke

CD163 Polyklonaler Antikörper

Katalog-Nr.: CL594-16646

1 Publikationen



Allgemeine Informationen

Katalog-Nr.: CL594-16646	GenBank-Zugangsnummer: BC051281	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 100ul, Konzentration: 1000 µg/ml von9332	GeneID (NCBI): 1000000000	Empfohlene Verdünnungen: IF 1:50-1:500
Nanodrop;	Vollständiger Name: CD163 molecule	Anregungs-/Emissionsmaxima- Wellenlängen: 588 nm / 604 nm
Wirt: Kaninchen	Berechnete Masse: 1156 aa, 125 kDa	
Isotyp: IgG		
Immunogen Katalognummer: AG10029		

Anwendungen

Geprüfte Anwendungen: IF	Positivkontrollen: IF : humanes Tonsillitisgewebe, humanes Lebergewebe
In Publikationen genannte Anwendungen: IF	
Getestete Reaktivität: Human	
Zitierte Arten: Human	

Hintergrundinformationen

CD163 is a transmembrane protein which belongs to the scavenger receptor cysteine-rich (SRCR) superfamily. This protein is a scavenger receptor for the hemoglobin-haptoglobin complex and is a marker for monocytes and macrophages. Soluble CD163 (sCD163), as a result of ectodomain shedding during inflammatory activation of macrophages, circulates in blood and has been suggested as a plasma/serum marker for macrophage activity.

Bemerkenswerte Veröffentlichungen

Verfasser	Pubmed ID	Journal	Anwendung
Yuyang Liu	35350788	Front Immunol	IF

Lagerung


Lagerungsbedingungen:
Bei -20°C lagern. Vor Licht schützen. Nach dem Versand ein Jahr stabil.
Lagerungspuffer:
BS mit 50% Glycerin, 0,05% Proclin300, 0,5% BSA, pH 7,3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Immunofluorescent analysis of (4% PFA) fixed human tonsillitis tissue using CoraLite®594 CD163 antibody (CL594-16646) at dilution of 1:200.



Immunofluorescent analysis of (4% PFA) fixed human tonsillitis tissue using CoraLite®594 CD163 antibody (CL594-16646) at dilution of 1:200.