

Nur für Forschungszwecke

CD86 Monoklonaler Antikörper

Katalog-Nr.:CL594-65068



Allgemeine Informationen

Katalog-Nr.:

CL594-65068

Größe:

100ug, 0.5 mg/ml

Wirt:

Ratte

Isotyp:

IgG2a, kappa

GenBank-Zugangsnummer:

BC013807

GeneID (NCBI):

12524

Vollständiger Name:

CD86 antigen

Reinigungsmethode:

Affinitätsreinigung

CloneNo.:

GL1

Anregungs-/Emissionsmaxima-
Wellenlängen:

588 nm / 604 nm

Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Maus

Hintergrundinformationen

CD86 (also known as B7-2) is a costimulatory molecule belonging to the immunoglobulin superfamily. Primarily expressed on antigen-presenting cells (APCs), including B cells, dendritic cells, and macrophages, CD86 is the ligand for two proteins at the cell surface of T cells, CD28 antigen and CTLA-4. Binding of CD86 with CD28 antigen is a costimulatory signal for activation of the T-cell. Binding of CD86 with CTLA-4 negatively regulates T-cell activation and diminishes the immune response. (PMID: 7513726; 1847722; 11029388; 27591335)

Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ untreated (dashed line) or LPS-treated mouse splenocytes were surface stained with 0.2 ug CoraLite® 594 Anti-Mouse CD86 (CL594-65068, Clone: GL1) (red), or 0.2 ug Isotype Control (blue). Cells were not fixed.



1X10⁶ untreated or LPS-treated mouse splenocytes were surface stained with 0.2 ug CoraLite® 594 Anti-Mouse CD86 (CL594-65068, Clone: GL1). Cells were not fixed.