

Nur für Forschungszwecke

CD94 Monoklonaler Antikörper

Katalog-Nr.:CL594-65240



Allgemeine Informationen

Katalog-Nr.: CL594-65240	GenBank-Zugangsnummer: BC120854	Reinigungsmethode: N/A
Größe: 100ug, 0.5 mg/ml	GeneID (NCBI): 16643	CloneNo.: 18D3
Wirt: Ratte	Vollständiger Name: killer cell lectin-like receptor, subfamily D, member 1	Anregungs-/Emissionsmaxima- Wellenlängen: 588 nm / 604 nm
Isotyp: IgG2a, kappa		

Anwendungen

Geprüfte Anwendungen:
FC

Getestete Reaktivität:
Maus

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:
PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ mouse splenocytes were surface stained with 0.5 ug CoraLite® 594 Anti-Mouse CD94 (CL594-65240, Clone: 18D3) or 0.5 ug Isotype Control. Cells were not fixed.