

Allgemeine Informationen

Katalog-Nr.:

CL594-66006

GenBank-Zugangsnummer:

GeneID (NCBI):

Größe:

100ul, Konzentration: 1000 µg/ml von
Nanodrop;

Vollständiger Name:

Wirt:

Maus

Isotyp:

IgG1

Reinigungsmethode:

Protein-G-Reinigung

CloneNo.:

1F5C6

Empfohlene Verdünnungen:

IF 1:50-1:500

Anregungs-/Emissionsmaxima-

Wellenlängen:

588 nm / 604 nm

Anwendungen

Geprüfte Anwendungen:

IF

Getestete Reaktivität:

Rekombinantes Protein

Positivkontrollen:

IF : Transfizierte HEK-293-Zellen,

Hintergrundinformationen

Protein tags are a protein or peptide sequences located either on the C- or N- terminal of the target protein, which can facilitate solubility, detection, purification, localization, and expression of the target protein. The HA tag corresponds to amino acid residues YPYDVPDYA of a surface glycoprotein -human influenza virus hemagglutinin (HA). The HA tag is commonly used for a variety of research applications including chromatin immunoprecipitation, ELISA, flow cytometry, western blotting, immunocytochemistry/immunofluorescence among of others.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Vor Licht schützen. Nach dem Versand ein Jahr stabil.

Lagerungspuffer:

BS mit 50% Glycerin, 0,05% Proclin300, 0,5% BSA, pH 7,3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Immunofluorescent analysis of (-20°C Ethanol)
fixed Transfected HEK-293 cells using
CoraLite®594 HA Tag antibody (CL594-66006,
Clone: 1F5C6) at dilution of 1:200.