

PTBP1 Monoklonaler Antikörper

Katalog-Nr.:CL594-67462

Allgemeine Informationen

Katalog-Nr.: CL594-67462	GenBank-Zugangsnummer: BC013694	Reinigungsmethode: Protein-G-Reinigung
Größe: 100ul, Konzentration: 1000 µg/ml von5725	GeneID (NCBI): von5725	CloneNo.: 3D6D8
Nanodrop;	Vollständiger Name: polypyrimidine tract binding protein	Empfohlene Verdünnungen: IF 1:50-1:500
Wirt: Maus	1	Anregungs-/Emissionsmaxima- Wellenlängen: 588 nm / 604 nm
Isotyp: IgG1	Berechnete Masse: 59 aa, 6 kDa	
Immunogen Katalognummer: AG28404	Beobachtete Masse: 60 kDa	

Anwendungen

Geprüfte Anwendungen: FC (Intra), IF	Positivkontrollen: IF : A375-Zellen,
Getestete Reaktivität: Human, Maus, Ratte	

Hintergrundinformationen


Polypyrimidine tract-binding protein (PTBP1), also known as heterogeneous nuclear ribonucleoprotein type I (hnRNP I), is a nuclear protein that binds pre-mRNAs in specific regions of the hnRNA-protein complexes sensitive to micrococcal nuclease. It involves in pre-mRNA splicing and the regulation of alternative splicing events. PTBP1 promotes RNA looping when bound to two separate polypyrimidine tracts in the same pre-mRNA, and promotes the binding of U2 snRNP to pre-mRNA. Besides, it can cooperate with RAVER1 to modulate switching between mutually exclusive exons during maturation of the TPM1 pre-mRNA, and represses the splicing of MAPT/Tau exon 10.

Lagerung


Lagerungsbedingungen:
Bei -20°C lagern. Vor Licht schützen. Nach dem Versand ein Jahr stabil.
Lagerungspuffer:
BS mit 50% Glycerin, 0,05% Proclin300, 0,5% BSA, pH 7,3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

Ausgewählte Validierungsdaten



Immunofluorescent analysis of (4% PFA) fixed A375 cells using CoraLite®594 PTBP1 antibody (CL594-67462, Clone: 3D6D8) at dilution of 1:200.



1X10⁶ HeLa cells were intracellularly stained with 0.4 ug CoraLite®594 Anti-Human PTBP1 (CL594-67462, Clone:3D6D8) (red), or 0.4 ug Control Antibody. Cells were fixed and permeabilized with Transcription Factor Staining Buffer Kit (PF00011).