

Nur für Forschungszwecke

CD48 Monoklonaler Antikörper

Katalog-Nr.:CL647-65227



Allgemeine Informationen

Katalog-Nr.: CL647-65227	GenBank-Zugangsnummer: BC060977	Reinigungsmethode: Protein-A-Reinigung
Größe: 100ug , 500 µg/ml	GeneID (NCBI): 12506	CloneNo.: HM48-1
Wirt: Armenian Hamster	Vollständiger Name: CD48 antigen	Anregungs-/Emissionsmaxima-Wellenlängen: 654 nm / 674 nm
Isotyp: IgG		

Anwendungen

Geprüfte Anwendungen:
FC
Getestete Reaktivität:
Maus

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Lagerungspuffer:
PBS mit 0,1% Natriumazid und 0,5% BSA, pH 7,3.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ mouse splenocytes were surface stained with 1 ug CoraLite® Plus 647 Anti-Mouse CD48 (CL647-65227, Clone: HM48-1) (red), or 1 ug Control Antibody. Cells were not fixed.