

Nur für Forschungszwecke

CD49b Monoklonaler Antikörper

Katalog-Nr.:CL647-65252



Allgemeine Informationen

Katalog-Nr.: CL647-65252	GenBank-Zugangsnummer: NM_002203	Reinigungsmethode: Affinitätsreinigung
Größe: 100tests , 5 µl/test	GeneID (NCBI): 3673	CloneNo.: AK7
Wirt: Maus	Vollständiger Name: integrin, alpha 2 (CD49B, alpha 2 subunit of VLA-2 receptor)	Anregungs-/Emissionsmaxima-Wellenlängen: 654 nm / 674 nm
Isotyp: IgG1, kappa	Berechnete Masse: 126 kDa	

Anwendungen

Geprüfte Anwendungen:
FC

Getestete Reaktivität:
Human

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:
PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ human PBMCs were surface stained with 5 ul CoraLite® Plus 647 Anti-Human CD49b (CL647-65252, Clone: AK7). Cells were not fixed. Monocytes were gated.