

KLRG1 Monoklonaler Antikörper

Katalog-Nr.: CL750-65084

Allgemeine Informationen

Katalog-Nr.: CL750-65084	GenBank-Zugangsnummer: BC101953	Reinigungsmethode: Affinitätsreinigung
Größe: 100ug, 0.5 mg/ml	GeneID (NCBI): 50928	CloneNo.: 2F1
Wirt: Golden Syrian Hamster	Vollständiger Name: killer cell lectin-like receptor subfamily G, member 1	Anregungs-/Emissionsmaxima- Wellenlängen: 755 nm / 780 nm
Isotyp: IgG		

Anwendungen

Geprüfte Anwendungen:
FC
Getestete Reaktivität:
Maus


Hintergrundinformationen

KLRG1 (killer cell lectin-like receptor subfamily G, member 1) is a C-type lectin inhibitory receptor that contains an immunoreceptor tyrosine-based inhibitory motif (ITIM) in its cytoplasmic domain. KLRG1 is expressed by subsets of NK and T cells, existing both as a monomer and as a disulfide-linked homodimer (PMID: 10741410; 19604491). It is considered to be a cell differentiation marker for NK and T cells and is strongly induced by viral and other infections (PMID: 20373518). Through interactions with members of the cadherin family, KLRG1 plays an inhibitory role on NK- and T-cell function (PMID: 19009530; 20373518).

Lagerung

Lagerungsbedingungen:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Lagerungspuffer:
PBS with 0.09% sodium azide.

Ausgewählte Validierungsdaten



1X10⁶ mouse splenocytes were surface co-stained with FITC Anti-Mouse NK1.1 (CD161) and 0.5 μ g CoraLite® Plus 750 Anti-Mouse KLRG1 (CL750-65084, Clone: 2F1) or 0.5 μ g Isotype Control. Cells were not fixed.