

Nur für Forschungszwecke

CD11b Monoklonaler Antikörper

Katalog-Nr.: CL750-65229



Allgemeine Informationen

Katalog-Nr.: CL750-65229	GenBank-Zugangsnummer: NM_012711	Reinigungsmethode: Affinitätsreinigung
Größe: 100ug, 0.5 mg/ml	GeneID (NCBI): 25021	CloneNo.: OX42
Wirt: Maus	Vollständiger Name: integrin alpha M	Anregungs-/Emissionsmaxima- Wellenlängen: 755 nm / 780 nm
Isotyp: IgG2a, kappa	Berechnete Masse: 127 kDa	

Anwendungen

Geprüfte Anwendungen:
FC
Getestete Reaktivität:
Ratte

Hintergrundinformationen

Integrins are cell adhesion receptors that are heterodimers composed of non-covalently associated α and β subunits (PMID: 9779984). CD11b, also known as Integrin alpha M or CR3A, belongs to the integrin alpha chain family. CD11b forms an α/β heterodimer with CD18 (integrin β 2). CD11b/CD18 is implicated in various adhesive interactions of monocytes, macrophages and granulocytes as well as in mediating the uptake of complement-coated particles and pathogens (PMID: 9558116; 20008295). CD11b/CD18 is a receptor for the complement protein fragment iC3b, and is also a receptor for fibrinogen, factor X and ICAM1 (PMID: 2971974; 15485828).


Lagerung

Lagerungsbedingungen:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Lagerungspuffer:
PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ Wistar rat splenocytes were surface stained with 0.5 ug CoraLite® Plus 750 Anti-Rat CD11b (CL750-65229, Clone: OX42) or 0.5 ug Control Antibody. Cells were not fixed.