

Nur für Forschungszwecke

CD94 Monoklonaler Antikörper

Katalog-Nr.: CL750-65240



Allgemeine Informationen

Katalog-Nr.:

CL750-65240

Größe:

100 µg, 500 µg/ml

Wirt:

Ratte

Isotyp:

IgG2a, kappa

GenBank-Zugangsnummer:

BC120854

GeneID (NCBI):

16643

Vollständiger Name:

killer cell lectin-like receptor,
subfamily D, member 1

Reinigungsmethode:

Affinitätsreinigung

CloneNo.:

18D3

Anregungs-/Emissionsmaxima-
Wellenlängen:

755 nm / 780 nm

Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Maus

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1x10⁶ mouse splenocytes were surface stained with FITC Anti-Mouse CD49b (HMA2) (FITC-65245, Clone: HMA2) and 0.5 ug CoraLite® Plus 750 Anti-Mouse CD94 (CL750-65240, Clone: 18D3) or CoraLite® Plus 750 Rat IgG2a isotype control. Cells were not fixed.