

Nur für Forschungszwecke

CD49b Monoklonaler Antikörper

Katalog-Nr.:CL750-65245



Allgemeine Informationen

Katalog-Nr.:

CL750-65245

Größe:

100µg, 500 µg/ml

Wirt:

Armenian Hamster

Isotyp:

IgG

GenBank-Zugangsnummer:

BC065139

GeneID (NCBI):

16398

Vollständiger Name:

integrin alpha 2

Reinigungsmethode:

Affinitätsreinigung

CloneNo.:

HMA2

Anregungs-/Emissionsmaxima-
Wellenlängen:

755 nm / 780 nm

Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Maus

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ mouse splenocytes were surface stained with 0.5 ug CoraLite® Plus 750 Anti-Mouse CD49b (CL750-65245, Clone: HMA2) or 0.5 ug Control Antibody. Cells were not fixed.