

Nur für Forschungszwecke

CD49b Monoklonaler Antikörper

Katalog-Nr.: CL750-65252



Allgemeine Informationen

Katalog-Nr.:

CL750-65252

Größe:

100 tests, 5 µl/test

Wirt:

Maus

Isotyp:

IgG1, kappa

GenBank-Zugangsnummer:

NM_002203

GeneID (NCBI):

3673

Vollständiger Name:

integrin, alpha 2 (CD49B, alpha 2 subunit of VLA-2 receptor)

Berechnete Masse:

126 kDa

Reinigungsmethode:

Affinitätsreinigung

CloneNo.:

AK7

Anregungs-/Emissionsmaxima-

Wellenlängen:

755 nm / 780 nm

Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Human

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:

PBS with 0.09% sodium azide and 0.5% BSA.

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ human PBMCs were surface stained with 5 ul CoraLite® Plus 750 Anti-Human CD49b (CL750-65252, Clone: AK7) or unstained. Cells were not fixed. Monocytes were gated.