

Nur für Forschungszwecke

CD18 Monoklonaler Antikörper

Katalog-Nr.:FITC-65178



Allgemeine Informationen

Katalog-Nr.:

FITC-65178

Größe:

100ug, 0.5 mg/ml

Wirt:

Ratte

Isotyp:

IgG2a, kappa

GenBank-Zugangsnummer:

BC145644

GeneID (NCBI):

16414

Vollständiger Name:

integrin beta 2

Reinigungsmethode:

Der gereinigte Antikörper ist mit Fluoresceinisothiocyanat (FITC) unter optimalen Bedingungen konjugiert. Das Konjugat wird mit Größenausschluss-Chromatographie gereinigt.

CloneNo.:

M18/2

Anregungs-/Emissionsmaxima-Wellenlängen:

494 nm / 520 nm

Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Maus

Hintergrundinformationen

CD18, also known as integrin beta-2, is a transmembrane protein that forms heterodimers with CD11a, CD11b, CD11c, CD11d (PMID: 1968349). CD11a/CD18 (LFA-1) is expressed by all leukocytes, while CD11b/CD18 (Mac-1) and CD11c/CD18 (p150,95) are normally restricted to expression on monocytes, macrophages, polymorphonuclear lymphocytes (PMNs), and natural killer cells (PMID: 1968349; 3109455; 10946284). CD18 plays an important role in mediating cell adhesion during immune response and defects of CD18 cause leukocyte adhesion deficiency (PMID: 3594570).

Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:

PBS mit 0,1% Natriumazid.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ mouse splenocytes were surface stained with 0.5 ug FITC Anti-Mouse CD18 (FITC-65178, Clone: M18/2) (red), or 0.5 ug Control Antibody. Cells were not fixed.



1X10⁶ mouse splenocytes were surface stained with 0.5 ug FITC Anti-Mouse CD18 (FITC-65178, Clone: M18/2) (red), or 0.5 ug Control Antibody. Cells were not fixed.