

Nur für Forschungszwecke

# STRO-1 Monoklonaler Antikörper

Katalog-Nr.: FITC-65184



## Allgemeine Informationen

Katalog-Nr.:

FITC-65184

Größe:

100 tests, 5 µl/test

Wirt:

Maus

Isotyp:

IgM, lambda

GenBank-Zugangsnummer:

GeneID (NCBI):

Vollständiger Name:

Reinigungsmethode:

Mit Präzipitation und  
Chromatographie aufgereinigt

CloneNo.:

STRO-1

Anregungs-/Emissionsmaxima-  
Wellenlängen:

495 nm / 524 nm

## Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Human

## Hintergrundinformationen

### Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:

PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Ausgewählte Validierungsdaten



1X10<sup>6</sup> human whole blood cells were surface stained with 5 ul FITC Plus Anti-Human STRO-1 (FITC-65184, Clone: STRO-1) or unstained. Cells were not fixed. CD45dim, CD34+ cells were gated.