

Nur für Forschungszwecke

CD48 Monoklonaler Antikörper

Katalog-Nr.:FITC-65227



Allgemeine Informationen

Katalog-Nr.:
FITC-65227

Größe:
100ug , 500 µg/ml

Wirt:
Armenian Hamster

Isotyp:
IgG

GenBank-Zugangsnummer:
BC060977

GeneID (NCBI):
12506

Vollständiger Name:
CD48 antigen

Reinigungsmethode:
Affinitätsreinigung

CloneNo.:
HM48-1

Anregungs-/Emissionsmaxima-
Wellenlängen:
495 nm / 524 nm

Anwendungen

Geprüfte Anwendungen:
FC

Getestete Reaktivität:
Maus

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ mouse splenocytes were surface stained with 1 ug FITC Plus Anti-Mouse CD48 (FITC-65227, Clone: HM48-1) (red), or 1 ug Control Antibody. Cells were not fixed.