

Nur für Forschungszwecke

Human IgG Kappa chain Monoklonaler Antikörper

Katalog-Nr.: HRP-67761



Allgemeine Informationen

Katalog-Nr.:

HRP-67761

GenBank-Zugangsnummer:

GeneID (NCBI):

Größe:

100ul, Konzentration: 1000 µg/ml von
Nanodrop;

Vollständiger Name:

Berechnete Masse:

155 kDa

Wirt:

Maus

Beobachtete Masse:

~25 kDa

Isotyp:

IgG1

Reinigungsmethode:

Protein-G-Reinigung

CloneNo.:

4D11F7

Empfohlene Verdünnungen:

WB 1:1000-1:4000

Anwendungen

Geprüfte Anwendungen:

WB

Getestete Reaktivität:

Human

Positivkontrollen:

WB : human plasma,

Hintergrundinformationen

This antibody detects the Kappa chain of human IgG,

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Vor Licht schützen.

Lagerungspuffer:

BS mit 50% Glycerin, 0,05% Proclin300, 0,5% BSA, pH 7,3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Ausgewählte Validierungsdaten



human plasma were subjected to SDS PAGE followed by western blot with HRP-67761 (Human IgG Kappa chain antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.