

Nur für Forschungszwecke

CD5 Monoklonaler Antikörper

Katalog-Nr.: PE-65152



Allgemeine Informationen

Katalog-Nr.: PE-65152	GenBank-Zugangsnummer: BC027901	Reinigungsmethode: Affinitätsreinigung
Größe: 100tests, 5 µl/test	GeneID (NCBI): 921	CloneNo.: UCHT2
Wirt: Maus	Vollständiger Name: CD5 molecule	Anregungs-/Emissionsmaxima-Wellenlängen: 496 nm, 565 nm / 578 nm
Isotyp: IgG1, kappa	Berechnete Masse: 495 aa, 55 kDa	

Anwendungen

Geprüfte Anwendungen:
FC

Getestete Reaktivität:
Human

Hintergrundinformationen

CD5 is a type I transmembrane glycoprotein of the scavenger receptor cysteine-rich family (PMID: 12403363). CD5 is expressed on a majority of thymocytes, mature T cells, B cell subsets, and peripheral blood dendritic cells (PMID: 9858516; 6156984; 10379049; 1384337). CD5 may act as a receptor in regulating T-cell proliferation. It functions as a negative regulator of TCR signaling during thymocyte development (PMID: 7542801).

Lagerung

Lagerungsbedingungen:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:
Phosphatbasierter Puffer mit 0,09% Natriumazid und 0,1% Gelatine, pH 7,2.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ human peripheral blood lymphocytes were surface stained with 0.5 ug PE-anti-human CD5 (PE-65152, clone UCHT2) (red) or 0.5 ug PE-mouse IgG1 isotype control (blue). Cells were not fixed.