

Nur für Forschungszwecke

CD137 Monoklonaler Antikörper

Katalog-Nr.: PE-65159



Allgemeine Informationen

Katalog-Nr.:

PE-65159

Größe:

100tests, 10 µl/test

Wirt:

Maus

Isotyp:

Mouse IgG1 kappa

GenBank-Zugangsnummer:

BC006196

GeneID (NCBI):

3604

Vollständiger Name:

tumor necrosis factor receptor
superfamily, member 9

Berechnete Masse:

28 kDa

Reinigungsmethode:

Der gereinigte Antikörper ist mit R-Phycoerythrin (PE) unter optimalen Bedingungen konjugiert. Das Konjugat wird mit Größenausschluss-Chromatographie gereinigt.

CloneNo.:

4B4-1

Anregungs-/Emissionsmaxima-Wellenlängen:

496 nm, 565 nm / 578 nm

Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Human, Nicht-Humane Primaten

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei 2-8°C lagern. Vor Licht schützen.

Lagerungspuffer:

PBS mit 0,1% Natriumazid.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ PHA-stimulated (5 µg/mL, 3 days) (red) or unstimulated (blue) human peripheral blood mononuclear cells (PBMCs) were surface stained with 10.00 µl PE Anti-Human CD137 (PE-65159, Clone: 4B4-1). Cells were not fixed.