

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-RBM14



Numéro de catalogue: 10196-1-AP

3 Publications

Informations de base

Numéro de catalogue:
10196-1-AP

Taille:
150ul, Concentration: 600 µg/ml by Nanodrop and 333 µg/ml by Bradford method using BSA as the standard;

Hôte:
Lapin

Isotype:
IgG

Immunogen Catalog Number:
AG0254

Numéro d'acquisition GenBank:
BC000488

Identification du gène (NCBI):
10432

Nom complet:
RNA binding motif protein 14

MW calculé
70 kDa

MW observés:
70-75 kDa

Méthode de purification:
Purification par affinité contre l'antigène

Dilutions recommandées:
WB 1:500-1:2000

Applications

Applications testées:
WB, ELISA

Demandes citées:
WB

Spécificité de l'espèce:
Humain

Espèces citées:
Humain

Contrôles positifs:

WB : cellules Jurkat, cellules HEK-293, cellules HeLa

Informations générales

RBM14 is a ribonucleoprotein that functions as a general nuclear coactivator, and an RNA splicing modulator. This protein contains two RNA recognition motifs (RRM) at the N-terminus, and multiple hexapeptide repeat domain at the C-terminus that interacts with thyroid hormone receptor-binding protein (TRBP), and is required for transcription activation.

Publications notables

Autrice	Pubmed ID	Journal	Application
Qianqian Yang	36434041	Exp Mol Med	WB
Yang Li	26979993	Brain Res	WB
Xinxin Zuo	33398173	Nat Struct Mol Biol	WB

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Jurkat cells were subjected to SDS PAGE followed by western blot with 10196-1-AP (RBM14 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.



Various lysates were subjected to SDS PAGE followed by western blot with 10196-1-AP (RBM14 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.