

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-Ajuba



Numéro de catalogue: 10653-1-AP

Phare

2 Publications

Informations de base

Numéro de catalogue:

10653-1-AP

Taille:

150ul, Concentration: 147 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG1011

Numéro d'acquisition GenBank:

BC007580

Identification du gène (NCBI):

84962

Nom complet:

jub, ajuba homolog (Xenopus laevis)

MW calculé

57 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Applications

Applications testées:

ELISA

Demandes citées:

IHC, WB

Spécificité de l'espèce:

Humain, rat, souris

Espèces citées:

Humain, souris

Informations générales

Ajuba belongs to the LIM domain proteins, which are involved in the assembly of the extracellular matrix and, along with associated proteins, regulate target genes that connect the extracellular matrix and the cytoskeleton. Ajuba consists of three LIM domains at the C-terminal and is located in the cytoplasm where it interacts with a variety of other proteins. Ajuba contributes to the formation or strengthening of cell-cell adhesion and regulates cell growth and differentiation decisions. Ajuba participates in the assembly of the extracellular matrix. Ajuba also plays an important role in some diseases.

Publications notables

Autrice	Pubmed ID	Journal	Application
Yandong Zhang	31731517	Molecules	WB
Ling Zhang	25839328	Am J Hum Genet	WB, IHC

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

