

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-Galectin 8



Numéro de catalogue: 10955-1-AP

Informations de base

Numéro de catalogue:	10955-1-AP	Numéro d'acquisition GenBank:	BC015818	Méthode de purification:	Purification par affinité contre l'antigène
Taille:	150ul, Concentration: 300 µg/ml by Nanodrop;	Identification du gène (NCBI):	3964	Dilutions recommandées:	WB 1:500-1:1000
Hôte:	Lapin	Nom complet:	lectin, galactoside-binding, soluble, 8		
Isotype:	IgG	MW calculé	34 kDa		
Immunogen Catalog Number:	AG1387	MW observés:	36 kDa, 40 kDa		

Applications

Applications testées:	WB, ELISA	Contrôles positifs:	WB : cellules HEK-293, cellules LNCaP
Spécificité de l'espèce:	Humain		

Informations générales

Galectin 8 is a member of the galectin family. Galectins are beta-galactoside-binding animal lectins with conserved carbohydrate recognition domains. The galectins have been implicated in many essential functions including development, differentiation, cell-cell adhesion, cell-matrix interaction, growth regulation, apoptosis, and RNA splicing. This gene is widely expressed in tumoral tissues and seems to be involved in integrin-like cell interactions

Stockage


Stockage:
Stocker à -20°C. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3
L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 10955-1-AP (Galectin 8 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.