

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-SF1



Numéro de catalogue: 11635-1-AP

1 Publications

Informations de base

Numéro de catalogue:

11635-1-AP

Taille:

150ul, Concentration: 900 µg/ml by Nanodrop;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG2209

Numéro d'acquisition GenBank:

BC020217

Identification du gène (NCBI):

7536

Nom complet:

splicing factor 1

MW calculé

638 aa, 68 kDa

MW observés:

65-71 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:500-1:3000

Applications

Applications testées:

WB, ELISA

Demandes citées:

WB

Spécificité de l'espèce:

Humain

Contrôles positifs:

WB : cellules HeLa, cellules HepG2

Informations générales

SF1 is a nuclear pre-mRNA splicing factor, it specifically recognizes the intron branch point sequence at the 3' splice site, together with the large subunit of U2 auxiliary factor (U2AF), and is required for the early stages of spliceosome assembly. It also plays a role in nuclear pre-mRNA retention and transcriptional repression. The molecular weight of SF1 ranges from 59 kDa to 71 kDa due to alternative splicing.

Publications notables

Autrice	Pubmed ID	Journal	Application
Rebecca Moschall	28166800	Retrovirology	WB
Morgan K Penny	37509222	Cancers (Basel)	IHC, IF

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 11635-1-AP (SF1 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.